eviFluor Duo Fluorometer

USER MANUAL

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Manufacturer:
Hombrechtikon Systems Engineering AG
Garstligweg 6
8634 Hombrechtikon
SWITZERLAND

This product may be covered by one or more patents
Document number: DMR-DOU-21200611-05122
Version: 002
Published on November 28, 2025

For more information, see www.on-deck-fluorometer.com

Table of Contents

USER MANUAL (ENGLISH)	1
1. Introduction	2
2. Conventions Used in This Manual	2
3. Safety Precautions and Hazards	3
4. Product Description	4
5. Getting Started	7
6. Maintenance and Cleaning	9
7. Service and Support	10
8. Self-Test and Status LED Indicators	11
9. Disposal	11
BENUTZERHANDBUCH (DEUTSCH)	13
MANUEL UTILISATEUR (FRANÇAIS)	25
用户手册 (简体中文)	37
ユーザーマニュアル (日本語)	47

1 Introduction

This user manual provides essential safety information, details on intended use, and installation instructions for the eviFluor Duo Fluorometer. This on-deck instrument functions exclusively when integrated into a liquid handling system. For information on installation and usage in conjunction with the liquid handling system, please refer to the documentation provided by the liquid handling system manufacturer.

Please note that the images provided in this manual may differ slightly from the actual instrument.

2 Conventions Used in This Manual

The following conventions are used throughout this manual:

WARNING

The term **WARNING** is used to indicate situations that could result in personal injury.

CAUTION

The term **CAUTION** is used to indicate situations that could result in damage to the instrument or to other equipment.

The information provided in this manual is intended to supplement, not replace, the standard safety requirements prevailing in the user's country.

3 Safety Precautions and Hazards

Read all product manuals and consult with trained personnel before operating the instrument. Always follow product labeling and manufacturer's recommendations. If in doubt as to how to proceed in any situation, contact your local liquid handler system provider or distributor.

**WARNING****CAUTION****Risk of personal injury and material damage**


Improper use of the eviFluor Duo Fluorometer may result in personal injury or damage to the instrument. The instrument must only be operated by qualified personnel with appropriate training. Always wear appropriate protective gear, including protective eyewear, gloves and suitable laboratory attire. Operate the instrument only when the liquid handling system safety shield is in place.

**WARNING****Risk of fire/explosion**

Do not allow cleaning fluid or decontamination agents to come into contact with the instrument's electrical parts. When using alcohol-based disinfectants for cleaning, ensure the instrument is switched off and allow flammable vapors time to disperse before turning it back on.

**CAUTION****Damage to the instrument**

Avoid spilling water or chemicals on the instrument.

4 Product Description

4.1 Scope of delivery

Item	Description
eviFluor Duo Fluorometer	Instrument
eviEtte Micro Cuvettes	Blister containing 5 eviEtte Micro Cuvette racks (total of 480 eviEtte Micro Cuvettes)
USB-C cable	3-m standard USB cable
User Manual	Multilingual User Manual
Cuvette Guide eviFluor Duo	Quantity: 2 pieces

4.2 Instrument description

No.	Name	Function
1	Status LED	Shows instrument status: Green: Idle status Red: Error Blue: Busy
2	eviEtte Micro Cuvettes	Transparent cuvettes
3	eviEtte Micro Cuvette rack	Rack containing 96 eviEtte Micro Cuvettes
4	Snapper	Snapper to release eviEtte Micro Cuvette racks from the eviFluor Duo Fluorometer.
5	Position of cuvette rack	Defines position of cuvette rack and cuvettes
6	Orientation of cuvette rack	Symbol of cuvette rack snapper indicating correct orientation to place cuvette rack on instrument
7	Cuvette Guide	Positions eviEtte Micro Cuvettes during measurements
8	USB-C port	For connection of the eviFluor Duo to the liquid handler

4.3 Intended use

The eviFluor Duo Fluorometer is an optical laboratory instrument designed to measure the fluorescence of biological or chemical analytes in aqueous samples. It is specifically designed for use with eviEtte Micro Cuvettes and an automated liquid handling system.

As an on-deck instrument, the eviFluor Duo Fluorometer functions exclusively when integrated into a liquid handling system and cannot be operated independently.

The instrument is intended for general laboratory use and research applications. It is not designed, nor classified, as an in-vitro diagnostic device (IVD) and must not be used for in-vitro diagnostic testing.

4.4 Technical specifications

Measurement principle	Fluorescence measurement
Light source	LED peak 470 and 625 nm
Excitation wavelength	470 nm / 625 nm
Emission wavelength	520 nm / 680 nm
Detectors	Si-Photodiode
Min. sample volume	1 µL
Total measurement volume	10 µL
Measurement acquisition time	1 s
Calibration type	2 point standard
Compatible cuvettes	eviEtte Micro Cuvettes
Dimensions	Base Socket (W x D): SBS format, 128 x 85 mm (5 x 3.4 in) H: 68 mm (2.7 in) including racks filled with cuvettes
Weight	0.6 kg (1.3 lb)
Communication interface	USB Type C (Protocol USB 2.0)
Environmental conditions ¹	15 - 30 °C (59 - 86 F) 20 - 70% relative humidity (non-condensing)

Altitude	Up to 2000 m (6500 ft)
Installation site	Indoor use only
Transport conditions	-10 - 60 °C (14 - 140 F)
Operating voltage	5 VDC (USB Type C)
Overvoltage category	I
Pollution degree	2
Certifications and compliance	CE marked (conforming with applicable EU directives) cMETus certified (conforming with US and Canadian safety standards)

¹ Refer to assay kit documentation for specific temperature requirements

5 Getting Started

For detailed instructions on installing and operating the eviFluor Duo Fluorometer on your liquid handling system, please refer to the documentation provided by the manufacturer of your liquid handling system.

The eviFluor Duo Fluorometer should only be mounted on a liquid handling system which is compatible with the eviFluor Duo Fluorometer. Avoid direct sunlight during operation. Ensure that all safety measures and procedures required by the liquid handling system are followed when operating the eviFluor Duo Fluorometer.

To maintain optimal performance:

- Wear gloves when handling the eviEtte Micro Cuvettes and racks to prevent accidental contamination of the eviEtte Micro Cuvettes.
- Handle eviEtte Micro Cuvette racks with care to prevent eviEtte Micro Cuvettes from falling out.
- Store eviEtte Micro Cuvettes protected from dust.
- Do not re-use eviEtte Micro Cuvettes.

5.1 Preparation of the system

Prepare the eviFluor Duo Fluorometer:

- Unpack the eviFluor Duo Fluorometer.
- Unpack the eviFluor Duo Cuvette Guide.
- Place the eviFluor Duo Cuvette Guide into the designated position on the eviFluor Duo Fluorometer. (Ensure that the nose on the cuvette guide fits into the slit on the instrument.)

Prepare the liquid handling system:

Connect the USB-A plug of the USB cable to the control unit of the liquid handling system and route the cable into the liquid handler.

5.2 Mounting an eviFluor Duo Fluorometer on the workdeck

- Choose an empty position on the workdeck of the liquid handling system to mount the eviFluor Duo Fluorometer.
- Route the USB cable to the selected position.
- Turn the eviFluor Duo Fluorometer upside down to connect the USB cable.
- Return the eviFluor Duo Fluorometer to its correct position.
- Place the eviFluor Duo Fluorometer on the designated position of the workdeck, ensuring it is correctly oriented as specified in the liquid handling system's software.
- Ensure that the USB cable does not interfere with the operation of the liquid handling system.

5.3 Preparation for operation

Prepare the eviFluor Duo Fluorometer for operation as follows:

- Remove empty racks by activating the snapper. Take the required number of eviEtte Micro Cuvette racks from their packaging and mount them onto the corresponding positions of the eviFluor Duo Fluorometer. Note the symbol engraved shows the correct orientation of the rack. (See Chapter 4.2 "Orientation of cuvette rack").

- Ensure that all eviEtte Micro Cuvettes are correctly positioned in the eviEtte Micro Cuvette racks.
- Ensure that the USB cable does not interfere with the operation of the liquid handling system.
- Once the instrument is powered on, an instrument self-test will automatically take place. Verify that the LED is green after the test.

5.4 Removal of eviFluor Duo Fluorometer from the workdeck

Remove the eviFluor Duo Fluorometer from the workdeck as follows:

- Remove any installed eviEtte Micro Cuvette racks to prevent eviEtte Micro Cuvettes from falling out. Store the cuvettes protected from dust. If applicable, store them protected from dust.
- Lift the eviFluor Duo Fluorometer and turn it upside down to disconnect the USB cable.
- Remove the eviFluor Duo Fluorometer from the liquid handling system.
- If applicable, remove the USB cable according to the instructions provided by the manufacturer of the liquid handling system.

6 Maintenance and Cleaning

6.1 Maintenance

The instrument is maintenance-free for standard usage conditions.

If replacement of the eviFluor Duo Cuvette Guide is required, follow these steps:

Firmly pull the eviFluor Duo Cuvette Guide out of the eviFluor Duo Fluorometer and install new cuvette guide.

6.2 Cleaning

Spills on the instrument may be biohazardous. Remove all spills immediately.

Repeated use of alcohol-based solutions on the surfaces of the instrument may damage the instrument's surface finish. Use alcohol-based disinfectants only for potentially biohazardous spills, but not for routine cleaning. For routine cleaning (e.g. dust removal), use deionized water only.

Particles can enter the optical path of the instrument. Use disposable lint-free cloth or paper towels to clean the instrument.

Cleaning liquids may damage the electronics of the instrument. Avoid the use of excessive liquids and do not clean the USB plug of the instrument.

To clean the instrument surface, follow the procedure below:

- Remove the instrument from the liquid handling system as described above.
- If applicable, remove the eviFluor Duo Cuvette Guide
- Moisten a disposable lint-free cloth or paper towel with disinfectant or deionized water.
- Clean the surface using the moistened cloth or paper towel.
- Repeat as needed.
- Reinstall the instrument and cuvette guide as described above.

7 Service and Support

For service and support, please contact the local support of your liquid handling system supplier.

8 Self-Test and Status LED Indicators

The eviFluor Duo Fluorometer performs a self-test when powered on.

The status LED indicates the result of the self-test as follows:

- Blue LED – Self-test in progress
- Green LED – Self-test successful
- Red LED – Self-test failed

If the status LED is red, turn off the power to the eviFluor Duo Fluorometer and verify that the cuvette guide is empty and correctly inserted. After checking the cuvette guide, turn the power back on. If the self-test fails again, contact the local support of your liquid handling system supplier.

9 Disposal

In conformance with the European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE), this instrument must not be disposed of in domestic waste. This also applies to countries outside of the EU, per their specific requirements.

Please dispose of this product in accordance with local regulations at the collecting point specified for electrical and electronic equipment.

Thank you for your contribution to environmental protection.

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herewith declares that the product

eviFluor Duo Fluorometer (Product No. EV-0020)

is in conformity with the following European Directives

Low Voltage Directive (LVD)	2014/35/EU
Electromagnetic Compatibility Directive (EMC)	2014/30/EU
Restriction of Hazardous Substances Directive (RoHS II) and Amendment	2011/65/EU 2015/863/EU

and the following harmonized standards

Product Safety:
EN 61010-1:2010 + A1:2019

EMC:
**EN 61326-1:2013, Class B, Basic
Electromagnetic Environment**

Other applicable standards and regulations

IEC 61010-1:2010, IEC 61010-1:2010/AMD1:2016 (IEC)
UL 61010-1:2012 +Am. 1: 2019 (US)
CSA C22.2 No 61010-1:2012 (CA)
JIS C 1010-1(2019) (JP)

IEC 61326-1:2020, Class B, Basic
Electromagnetic Environment (IEC)
EN IEC 61326-1:2021, Class B, Basic
Electromagnetic Environment (CENELEC)
47 CFR Part 15 Subpart B, Class B (US)
ICES-001 Issue 5, Class B (CA)

Hombrechtikon, 13. May 2025

Michael Collasius
CEO